

## **Tetradecanoic acid-ethyl ester**

Art. ID                      LM24-OP-4008010  
Unit                         1 mL  
Deliverydetails            No Dangerous Good /not restricted

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter              | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------|------------|---------------------|------|--------|--------|
|                  | Tetradecanoic acid-ethyl ester | [124-06-1] |                     |      |        |        |