

Docosanoic acid-methyl ester

Art. ID LM24-OP-4009276
Unit 100 mg
Deliverydetails No Dangerous Good /not restricted

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|------------|---------------------|------|--------|--------|
| | Docosanoic acid methyl ester (Behenic acid methyl ester) | [929-77-1] | | | | |