

2-sec-Butyl-3-methoxypyrazine

Art. ID LM24-OP-4009435
Unit 100 mg
Deliverydetails De Minimis UN3334

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------------|--------------|---------------------|------|--------|--------|
| | 2-sec-Butyl-3-methoxy pyrazine | [24168-70-5] | | | | |