

**Bis-(4-aminophenyl)-methane (4,4'-Methylenedianiline)**

Art. ID LM24-OP-4009699  
Unit 250 mg  
Deliverydetails De Minimis UN2651 III E1

Description

Shelf-life: 24 months

| Text/Information | Analyte/Parameter   | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|------------|---------------------|------|--------|--------|
|                  | Bis-(4-aminophenyl)-met<br>hane (4,4'-Methylenedia<br>niline) | [101-77-9] |                     |      |        |        |