

Chromium (III) 1000 mg/l in diluted Nitric Acid

Art. ID LM24-OP-4010002
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | Chromium (III) | | 1000 | mg/l | | |