

**Iron in Flat Glass, slide approx. 100 mm x 50 mm, 3,2 mm thickness**

Art. ID	BAM-S050
Unit	each
Deliverydetails	No Dangerous Good /not restricted

**Description**

The CRM is intended for quality control. The minimum sample size for wet chemical analysis is 0.2 g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified value	Fe(II)		0,0027 ± 0,0005	%		
Certified value	Fe(total)		0,0084 ± 0,0012	%		
Values for information	Fe(III), calculated		0,0057 ± 0,0013	%		
Values for information	Fe(II), calculated as F e2O3		0,0038 ± 0,0008	%		