

Iron in Flat Glass, slide approx. 100 mm x 50 mm, 3,2 mm thickness

| | |
|-----------------|-----------------------------------|
| Art. ID | BAM-S050 |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

The CRM is intended for quality control. The minimum sample size for wet chemical analysis is 0.2 g

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------------|---------------------------------|------------|---------------------|------|--------|--------|
| Certified value | Fe(II) | | 0,0027 ± 0,0005 | % | | |
| Certified value | Fe(total) | | 0,0084 ± 0,0012 | % | | |
| Values for information | Fe(III), calculated | | 0,0057 ± 0,0013 | % | | |
| Values for information | Fe(II), calculated as F e2O3 | | 0,0038 ± 0,0008 | % | | |