

**Iron in Flat Glass, slide approx. 100 mm x 50 mm, 5,9 mm thickness**

Art. ID	BAM-S051
Unit	each
Deliverydetails	No Dangerous Good /not restricted

**Description**

The CRM is intended for quality control. The minimum sample size for wet chemical analysis is 0.2 g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified value	Fe(II)		0.0158 ± 0.0015	%		
Certified value	Fe(total)		0.0481 ± 0.0017	%		
Values for information	Fe(III), calculated		0.0323 ± 0.0022	%		
Values for information	Fe(II), calculated as F e2O3		0.0226 ± 0.0022	%		