

**Iron in Flat Glass, slide approx. 100 mm x 50 mm, 3,9 mm thickness**

Art. ID	BAM-S052
Unit	each
Deliverydetails	No Dangerous Good /not restricted

**Description**

The CRM is intended for quality control. The minimum sample size for wet chemical analysis is 0.2 g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified value	Fe(II)		0.164 ± 0.009	%		
Certified value	Fe(total)		0.597 ± 0.011	%		
Values for information	Fe(III), calculated		0.433 ± 0.014	%		
Values for information	Fe(II), calculated as F e2O3		0.234 ± 0.013	%		