

Cylindrospermopsin standard solution in deionized water

Art. ID NRC-CRM-CYN
Unit 0,5 mL
Deliverydetails No Dangerous Good (Dry ice shipment) /not restricted

Description

Each ampoule contains approximately 0.5 mL of solution with 30 ± 2 µM of cylindrospermopsin in filtered, deionized water

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cylindrospermopsin	[143545-90-8]	30 ± 2	µM		