

Cylindrospermopsin standard solution in deionized water

Art. ID NRC-CRM-CYN
Unit 0,5 mL
Deliverydetails Dry ice shipment /not restricted

Description

Each ampoule contains approximately 0.5 mL of solution with $30 \pm 2 \mu\text{M}$ of cylindrospermopsin in filtered, deionized water

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cylindrospermopsin	[143545-90-8]	30 ± 2	μM		