

### **Cylindrospermopsin standard solution in deionized water**

Art. ID                      NRC-CRM-CYN  
Unit                         0,5 mL  
Deliverydetails            Dry ice shipment /not restricted

#### **Description**

Each ampoule contains approximately 0.5 mL of solution with  $30 \pm 2 \mu\text{M}$  of cylindrospermopsin in filtered, deionized water

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cylindrospermopsin	[143545-90-8]	$30 \pm 2$	$\mu\text{M}$		