

## **Stainless Steel, PH13-8Mo/UNS S13800, chips**

Art. ID IARM-21D-C  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1,03	%		
	Carbon (C)	[7440-44-0]	0,032	%		
	Cobalt (Co)	[7440-48-4]	0,078	%		
	Chromium (Cr)	[7440-47-3]	12,69	%		
	Copper (Cu)	[7440-50-8]	0,017	%		
	Manganese (Mn)	[7439-96-5]	0,052	%		
	Molybdenum (Mo)	[7439-98-7]	2,23	%		
	Nitrogen (N)	[7727-37-9]	0,0037	%		
	Niobium (Nb)	[7440-03-1]	~0,005	%		
	Nickel (Ni)	[7440-02-0]	8,29	%		
	Phosphorus (P)	[7723-14-0]	0,008	%		
	Sulfur (S)	[7704-34-9]	~0,0014	%		
	Silicon (Si)	[7440-21-3]	0,039	%		
	Titanium (Ti)	[7440-32-6]	0,016	%		
	Vanadium (V)	[7440-62-2]	0,017	%		
	Tungsten (W)	[7440-33-7]	~0,012	%		