

## Titanium, Ti 6-7, chips

Art. ID IARM-300A-C  
Unit 65 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	5,97	%		
	Carbon (C)	[7440-44-0]	0,006	%		
	Chromium (Cr)	[7440-47-3]	0,009	%		
	Copper (Cu)	[7440-50-8]	~0,003	%		
	Iron (Fe)	[7439-89-6]	0,19	%		
	Hydrogen (H)	[1333-74-0]	0,0019	%		
	Manganese (Mn)	[7439-96-5]	~0,001	%		
	Molybdenum (Mo)	[7439-98-7]	~0,001	%		
	Nitrogen (N)	[7727-37-9]	0,004	%		
	Niobium (Nb)	[7440-03-1]	6,9	%		
	Nickel (Ni)	[7440-02-0]	0,008	%		
	Oxygen (O)	[7782-44-7]	0,16	%		
	Silicon (Si)	[7440-21-3]	0,019	%		
	Tin (Sn)	[7440-31-5]	~0,01	%		
	Titanium (Ti)	[7440-32-6]	~86,7	%		
	Vanadium (V)	[7440-62-2]	~0,01	%		