

Stainless & High Temperature Steel, Alloy C-300/UNS K93120, disc 38 x 3 mm

Art. ID IARM-99D-X
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,117	%		
	Boron (B)	[7440-42-8]	0,0026	%		
	Carbon (C)	[7440-44-0]	~0,006	%		
	Cobalt (Co)	[7440-48-4]	9,24	%		
	Chromium (Cr)	[7440-47-3]	~0,12	%		
	Copper (Cu)	[7440-50-8]	~0,045	%		
	Manganese (Mn)	[7439-96-5]	~0,013	%		
	Molybdenum (Mo)	[7439-98-7]	4,8	%		
	Nitrogen (N)	[7727-37-9]	~0,0014	%		
	Niobium (Nb)	[7440-03-1]	~0,011	%		
	Nickel (Ni)	[7440-02-0]	18,4	%		
	Phosphorus (P)	[7723-14-0]	~0,004	%		
	Sulfur (S)	[7704-34-9]	0,0011	%		
	Silicon (Si)	[7440-21-3]	~0,03	%		
	Titanium (Ti)	[7440-32-6]	0,67	%		
	Vanadium (V)	[7440-62-2]	~0,037	%		
	Tungsten (W)	[7440-33-7]	~0,010	%		