

Copper base, Type: Brass - Silicon, CRM, CDA 875, Chips 100 g

Art. ID IARM-151B-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,002	%		
	Carbon (C)	[7440-44-0]	0,005	%		
	Copper (Cu)	[7440-50-8]	84	%		
	Iron (Fe)	[7439-89-6]	0,025	%		
	Manganese (Mn)	[7439-96-5]	0,002	%		
	Nickel (Ni)	[7440-02-0]	0,011	%		
	Phosphorus (P)	[7723-14-0]	0,003	%		
	Lead (Pb)	[7439-92-1]	0,013	%		
	Sulfur (S)	[7704-34-9]	<0,001	%		
	Silicon (Si)	[7440-21-3]	3,11	%		
	Tin (Sn)	[7440-31-5]	0,009	%		
	Zinc (Zn)	[7440-66-6]	12,94	%		