

**Copper base, Type: Silver Alloy, RM, NARloy-Z, Chips 100 g**

Art. ID IARM-160A-C  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	3,03	%		
	Aluminium (Al)	[7429-90-5]	(<0,01)	%		
	Carbon (C)	[7440-44-0]	0,003	%		
	Cobalt (Co)	[7440-48-4]	(<0,01)	%		
	Chromium (Cr)	[7440-47-3]	(<0,01)	%		
	Iron (Fe)	[7439-89-6]	(<0,01)	%		
	Manganese (Mn)	[7439-96-5]	(<0,01)	%		
	Nickel (Ni)	[7440-02-0]	(<0,01)	%		
	Phosphorus (P)	[7723-14-0]	~0,004	%		
	Lead (Pb)	[7439-92-1]	(<0,01)	%		
	Sulfur (S)	[7704-34-9]	(<0,003)	%		
	Silicon (Si)	[7440-21-3]	(<0,01)	%		
	Tin (Sn)	[7440-31-5]	(<0,01)	%		
	Zinc (Zn)	[7440-66-6]	(<0,01)	%		
	Zirconium (Zr)	[7440-67-7]	0,4	%		