

Copper base, Type: Bronze - Tin, CRM, 941, Disc 31 mm x 18 mm

Art. ID IARM-184A-G
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0016	%		
	Arsenic (As)	[7440-38-2]	0,01	%		
	Nickel (Ni)	[7440-02-0]	0,3	%		
	Phosphorus (P)	[7723-14-0]	0,008	%		
	Lead (Pb)	[7439-92-1]	19	%		
	Sulfur (S)	[7704-34-9]	0,021	%		
	Antimony (Sb)	[7440-36-0]	0,27	%		
	Tin (Sn)	[7440-31-5]	6	%		
	Zinc (Zn)	[7440-66-6]	0,37	%		