

Nickel Alloy, RM 200, HP Ni-High, RoHS Monitor, Chips 100 g

Art. ID IARM-190A-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,00109	%		
	Aluminium (Al)	[7429-90-5]	0,005	%		
	Arsenic (As)	[7440-38-2]	0,0028	%		
	Boron (B)	[7440-42-8]	<0,0005	%		
	Beryllium (Be)	[7440-41-7]	(<0,0001)	%		
	Bismuth (Bi)	[7440-69-9]	0,00111	%		
	Carbon (C)	[7440-44-0]	0,0022	%		
	Calcium (Ca)	[7440-70-2]	<0,0010	%		
	Cadmium (Cd)	[7440-43-9]	0,0005	%		
	Cobalt (Co)	[7440-48-4]	0,0008	%		
	Chromium (Cr)	[7440-47-3]	~0,0001	%		
	Copper (Cu)	[7440-50-8]	0,0017	%		
	Iron (Fe)	[7439-89-6]	0,0099	%		
	Gallium (Ga)	[7440-55-3]	<0,00005	%		
	Magnesium (Mg)	[7439-95-4]	~0,0006	%		
	Manganese (Mn)	[7439-96-5]	0,00018	%		
	Molybdenum (Mo)	[7439-98-7]	(<0,0001)	%		
	Nitrogen (N)	[7727-37-9]	~0,0001	%		
	Nickel (Ni)	[7440-02-0]	rem	%		
	Oxygen (O)	[7782-44-7]	~0,0019	%		
	Phosphorus (P)	[7723-14-0]	0,0034	%		
	Lead (Pb)	[7439-92-1]	0,00093	%		
	Sulfur (S)	[7704-34-9]	0,00033	%		
	Antimony (Sb)	[7440-36-0]	0,0011	%		
	Selenium (Se)	[7782-49-2]	0,00065	%		
	Silicon (Si)	[7440-21-3]	0,0028	%		
	Tin (Sn)	[7440-31-5]	0,00062	%		
	Tellurium (Te)	[13494-80-9]	0,00089	%		
	Titanium (Ti)	[7440-32-6]	~0,0006	%		
	Thallium (Tl)	[7440-28-0]	0,00058	%		
	Vanadium (V)	[7440-62-2]	(<0,00005)	%		

Zinc (Zn)

[7440-66-6]

0,00081

%