

Iron base, Steel - Low Alloy, CRM, 8740, Chips 100 g

Art. ID IARM-252E-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,028	%		
	Arsenic (As)	[7440-38-2]	0,0046	%		
	Carbon (C)	[7440-44-0]	0,413	%		
	Cobalt (Co)	[7440-48-4]	0,0093	%		
	Chromium (Cr)	[7440-47-3]	0,486	%		
	Copper (Cu)	[7440-50-8]	0,164	%		
	Manganese (Mn)	[7439-96-5]	0,87	%		
	Molybdenum (Mo)	[7439-98-7]	0,204	%		
	Nitrogen (N)	[7727-37-9]	0,0064	%		
	Nickel (Ni)	[7440-02-0]	0,407	%		
	Phosphorus (P)	[7723-14-0]	~0,009	%		
	Sulfur (S)	[7704-34-9]	~0,012	%		
	Silicon (Si)	[7440-21-3]	0,257	%		
	Tin (Sn)	[7440-31-5]	0,0075	%		
	Titanium (Ti)	[7440-32-6]	0,001	%		
	Vanadium (V)	[7440-62-2]	~0,0028	%		