

Copper base, Type: Sebiloy, CRM, Federalloy I-836, Chips 100 g

Art. ID IARM-265A-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	~0,002	%		
	Aluminium (Al)	[7429-90-5]	0,003	%		
	Arsenic (As)	[7440-38-2]	~0,005	%		
	Bismuth (Bi)	[7440-69-9]	2,4	%		
	Copper (Cu)	[7440-50-8]	~90	%		
	Iron (Fe)	[7439-89-6]	0,013	%		
	Nickel (Ni)	[7440-02-0]	0,69	%		
	Phosphorus (P)	[7723-14-0]	0,024	%		
	Lead (Pb)	[7439-92-1]	0,011	%		
	Antimony (Sb)	[7440-36-0]	0,015	%		
	Silicon (Si)	[7440-21-3]	0,003	%		
	Tin (Sn)	[7440-31-5]	4,4	%		
	Zinc (Zn)	[7440-66-6]	2,45	%		