

Copper base, Type: Bronze - Tin, CRM, 947, Disc 31 mm x 18 mm

Art. ID IARM-267A-G
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	~0,002	%		
	Aluminium (Al)	[7429-90-5]	0,003	%		
	Arsenic (As)	[7440-38-2]	~0,004	%		
	Bismuth (Bi)	[7440-69-9]	~0,005	%		
	Carbon (C)	[7440-44-0]	~0,003	%		
	Copper (Cu)	[7440-50-8]	87,8	%		
	Iron (Fe)	[7439-89-6]	0,019	%		
	Nickel (Ni)	[7440-02-0]	5,1	%		
	Phosphorus (P)	[7723-14-0]	0,037	%		
	Lead (Pb)	[7439-92-1]	0,026	%		
	Sulfur (S)	[7704-34-9]	0,0014	%		
	Antimony (Sb)	[7440-36-0]	<0,03	%		
	Selenium (Se)	[7782-49-2]	~0,002	%		
	Silicon (Si)	[7440-21-3]	0,003	%		
	Tin (Sn)	[7440-31-5]	4,95	%		
	Zinc (Zn)	[7440-66-6]	2,06	%		