

Iron base, Steel - High Alloy, CRM, A-286, Chips 100 g

Art. ID IARM-26D-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,29	%		
	Boron (B)	[7440-42-8]	0,0063	%		
	Carbon (C)	[7440-44-0]	0,038	%		
	Cobalt (Co)	[7440-48-4]	0,04	%		
	Chromium (Cr)	[7440-47-3]	14,29	%		
	Copper (Cu)	[7440-50-8]	0,047	%		
	Manganese (Mn)	[7439-96-5]	0,224	%		
	Molybdenum (Mo)	[7439-98-7]	1,23	%		
	Nitrogen (N)	[7727-37-9]	0,0035	%		
	Niobium (Nb)	[7440-03-1]	~0,007	%		
	Nickel (Ni)	[7440-02-0]	24,6	%		
	Phosphorus (P)	[7723-14-0]	0,013	%		
	Sulfur (S)	[7704-34-9]	~0,0008	%		
	Silicon (Si)	[7440-21-3]	~0,05	%		
	Tin (Sn)	[7440-31-5]	0,0039	%		
	Titanium (Ti)	[7440-32-6]	2,17	%		
	Vanadium (V)	[7440-62-2]	0,223	%		
	Tungsten (W)	[7440-33-7]	0,036	%		