

## **Iron base, Steel - Tool, CRM, T-4, Chips 100 g**

Art. ID IARM-281A-C  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,007	%		
	Arsenic (As)	[7440-38-2]	~0,02	%		
	Boron (B)	[7440-42-8]	~0,003	%		
	Carbon (C)	[7440-44-0]	0,74	%		
	Cobalt (Co)	[7440-48-4]	4,8	%		
	Chromium (Cr)	[7440-47-3]	3,89	%		
	Copper (Cu)	[7440-50-8]	0,096	%		
	Magnesium (Mg)	[7439-95-4]	~0,003	%		
	Manganese (Mn)	[7439-96-5]	0,3	%		
	Molybdenum (Mo)	[7439-98-7]	0,49	%		
	Nitrogen (N)	[7727-37-9]	0,0064	%		
	Niobium (Nb)	[7440-03-1]	0,094	%		
	Nickel (Ni)	[7440-02-0]	0,15	%		
	Oxygen (O)	[7782-44-7]	~0,003	%		
	Phosphorus (P)	[7723-14-0]	0,015	%		
	Sulfur (S)	[7704-34-9]	0,019	%		
	Silicon (Si)	[7440-21-3]	0,29	%		
	Tin (Sn)	[7440-31-5]	0,02	%		
	Titanium (Ti)	[7440-32-6]	0,004	%		
	Vanadium (V)	[7440-62-2]	0,9	%		
	Tungsten (W)	[7440-33-7]	17,6	%		
	Zirconium (Zr)	[7440-67-7]	~0,002	%		