

## **Iron Base, Steel - Maraging, CRM, C-350, Chips 100 g**

Art. ID IARM-309A-C  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,11	%		
	Boron (B)	[7440-42-8]	0,0032	%		
	Carbon (C)	[7440-44-0]	0,0059	%		
	Cobalt (Co)	[7440-48-4]	12,3	%		
	Chromium (Cr)	[7440-47-3]	0,053	%		
	Copper (Cu)	[7440-50-8]	0,023	%		
	Manganese (Mn)	[7439-96-5]	0,018	%		
	Molybdenum (Mo)	[7439-98-7]	4,71	%		
	Nitrogen (N)	[7727-37-9]	0,001	%		
	Niobium (Nb)	[7440-03-1]	0,004	%		
	Nickel (Ni)	[7440-02-0]	18,4	%		
	Oxygen (O)	[7782-44-7]	0,0005	%		
	Phosphorus (P)	[7723-14-0]	0,004	%		
	Sulfur (S)	[7704-34-9]	0,0006	%		
	Silicon (Si)	[7440-21-3]	0,02	%		
	Tin (Sn)	[7440-31-5]	~0,001	%		
	Tantalum (Ta)	[7440-25-7]	~0,006	%		
	Titanium (Ti)	[7440-32-6]	1,47	%		
	Vanadium (V)	[7440-62-2]	0,01	%		
	Tungsten (W)	[7440-33-7]	0,01	%		
	Zirconium (Zr)	[7440-67-7]	0,008	%		