

Nickel base, Haynes 230, CRM, Chips 100 g

Art. ID IARM-68E-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,3	%		
	Boron (B)	[7440-42-8]	0,007	%		
	Carbon (C)	[7440-44-0]	0,099	%		
	Cobalt (Co)	[7440-48-4]	0,16	%		
	Chromium (Cr)	[7440-47-3]	21,88	%		
	Copper (Cu)	[7440-50-8]	0,022	%		
	Iron (Fe)	[7439-89-6]	1,06	%		
	Lanthanum (La)	[7439-91-0]	~0,008	%		
	Magnesium (Mg)	[7439-95-4]	~0,006	%		
	Manganese (Mn)	[7439-96-5]	0,51	%		
	Molybdenum (Mo)	[7439-98-7]	1,18	%		
	Nitrogen (N)	[7727-37-9]	0,05	%		
	Niobium (Nb)	[7440-03-1]	0,035	%		
	Nickel (Ni)	[7440-02-0]	~59,9	%		
	Oxygen (O)	[7782-44-7]	0,0007	%		
	Phosphorus (P)	[7723-14-0]	~0,005	%		
	Sulfur (S)	[7704-34-9]	~0,0005	%		
	Silicon (Si)	[7440-21-3]	0,39	%		
	Titanium (Ti)	[7440-32-6]	0,015	%		
	Vanadium (V)	[7440-62-2]	~0,007	%		
	Tungsten (W)	[7440-33-7]	14,6	%		