

*****sold out*** Copper base, Copper, CRM, 110, Chips 100 g**

Art. ID IARM-70C-C
Unit 100 g (chips)
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	~0,0011	%		
	Aluminium (Al)	[7429-90-5]	~0,0014	%		
	Arsenic (As)	[7440-38-2]	~0,0009	%		
	Boron (B)	[7440-42-8]	~0,0004	%		
	Beryllium (Be)	[7440-41-7]	<0,002	%		
	Bismuth (Bi)	[7440-69-9]	~0,002	%		
	Carbon (C)	[7440-44-0]	~0,002	%		
	Cadmium (Cd)	[7440-43-9]	~0,0008	%		
	Cobalt (Co)	[7440-48-4]	~0,0014	%		
	Chromium (Cr)	[7440-47-3]	~0,0002	%		
	Copper (Cu)	[7440-50-8]	99,94	%		
	Iron (Fe)	[7439-89-6]	~0,0016	%		
	Germanium (Ge)	[7440-56-4]	<0,0002	%		
	Indium (In)	[7440-74-6]	<0,0001	%		
	Magnesium (Mg)	[7439-95-4]	~0,0003	%		
	Manganese (Mn)	[7439-96-5]	~0,0002	%		
	Nitrogen (N)	[7727-37-9]	<0,008	%		
	Niobium (Nb)	[7440-03-1]	~0,004	%		
	Nickel (Ni)	[7440-02-0]	~0,0004	%		
	Oxygen (O)	[7782-44-7]	~0,002	%		
	Phosphorus (P)	[7723-14-0]	~0,0014	%		
	Lead (Pb)	[7439-92-1]	~0,0013	%		
	Sulfur (S)	[7704-34-9]	0,0008	%		
	Antimony (Sb)	[7440-36-0]	~0,003	%		
	Selenium (Se)	[7782-49-2]	~0,001	%		
	Silicon (Si)	[7440-21-3]	~0,0006	%		
	Tin (Sn)	[7440-31-5]	0,0005	%		
	Tellurium (Te)	[13494-80-9]	~0,001	%		
	Zinc (Zn)	[7440-66-6]	~0,002	%		
	Zirconium (Zr)	[7440-67-7]	~0,001	%		