

**Copper base, Type: Brass - Leaded, CRM, 314, Disc 31 mm x 18 mm**

Art. ID IARM-72B-G  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	~0,002	%		
	Arsenic (As)	[7440-38-2]	~0,003	%		
	Carbon (C)	[7440-44-0]	0,002	%		
	Copper (Cu)	[7440-50-8]	90,08	%		
	Iron (Fe)	[7439-89-6]	0,007	%		
	Nitrogen (N)	[7727-37-9]	~0,0002	%		
	Nickel (Ni)	[7440-02-0]	0,004	%		
	Oxygen (O)	[7782-44-7]	0,002	%		
	Phosphorus (P)	[7723-14-0]	0,005	%		
	Lead (Pb)	[7439-92-1]	1,99	%		
	Sulfur (S)	[7704-34-9]	0,0015	%		
	Antimony (Sb)	[7440-36-0]	0,006	%		
	Silicon (Si)	[7440-21-3]	~0,002	%		
	Tin (Sn)	[7440-31-5]	0,029	%		
	Zinc (Zn)	[7440-66-6]	7,81	%		