

Copper base, Brass - Naval, CRM, 485, Chips 100 g

Art. ID IARM-76D-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0014	%		
	Aluminium (Al)	[7429-90-5]	~0,002	%		
	Arsenic (As)	[7440-38-2]	~0,004	%		
	Beryllium (Be)	[7440-41-7]	~0,0006	%		
	Bismuth (Bi)	[7440-69-9]	0,0011	%		
	Carbon (C)	[7440-44-0]	~0,002	%		
	Cadmium (Cd)	[7440-43-9]	0,001	%		
	Cobalt (Co)	[7440-48-4]	0,001	%		
	Chromium (Cr)	[7440-47-3]	0,0006	%		
	Copper (Cu)	[7440-50-8]	60,7	%		
	Iron (Fe)	[7439-89-6]	0,013	%		
	Magnesium (Mg)	[7439-95-4]	~0,0005	%		
	Manganese (Mn)	[7439-96-5]	0,0006	%		
	Nickel (Ni)	[7440-02-0]	~0,003	%		
	Oxygen (O)	[7782-44-7]	~0,0006	%		
	Phosphorus (P)	[7723-14-0]	0,0018	%		
	Lead (Pb)	[7439-92-1]	1,69	%		
	Sulfur (S)	[7704-34-9]	0,0012	%		
	Antimony (Sb)	[7440-36-0]	0,004	%		
	Silicon (Si)	[7440-21-3]	0,0037	%		
	Tin (Sn)	[7440-31-5]	0,73	%		
	Zinc (Zn)	[7440-66-6]	36,8	%		