

**Copper base, Type: Bronze - Phosphorus, CRM, 544, Disc 38 mm x 3 mm**

Art. ID IARM-78B-X  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,008	%		
	Aluminium (Al)	[7429-90-5]	~0,002	%		
	Arsenic (As)	[7440-38-2]	<0,003	%		
	Bismuth (Bi)	[7440-69-9]	~0,001	%		
	Carbon (C)	[7440-44-0]	<0,002	%		
	Cobalt (Co)	[7440-48-4]	<0,0005	%		
	Copper (Cu)	[7440-50-8]	87,7	%		
	Iron (Fe)	[7439-89-6]	0,02	%		
	Manganese (Mn)	[7439-96-5]	~0,002	%		
	Nickel (Ni)	[7440-02-0]	0,077	%		
	Phosphorus (P)	[7723-14-0]	0,19	%		
	Lead (Pb)	[7439-92-1]	3,87	%		
	Sulfur (S)	[7704-34-9]	0,01	%		
	Antimony (Sb)	[7440-36-0]	0,01	%		
	Silicon (Si)	[7440-21-3]	<0,002	%		
	Tin (Sn)	[7440-31-5]	4,73	%		
	Zinc (Zn)	[7440-66-6]	3,55	%		