

**Copper base, Type: Bronze - Aluminum, CRM, CDA 623, Chips 100 g**

Art. ID IARM-79C-C  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	<0,005	%		
	Aluminium (Al)	[7429-90-5]	9,2	%		
	Arsenic (As)	[7440-38-2]	~0,003	%		
	Carbon (C)	[7440-44-0]	~0,003	%		
	Cobalt (Co)	[7440-48-4]	<0,005	%		
	Chromium (Cr)	[7440-47-3]	~0,002	%		
	Copper (Cu)	[7440-50-8]	87,6	%		
	Iron (Fe)	[7439-89-6]	2,28	%		
	Manganese (Mn)	[7439-96-5]	0,2	%		
	Nitrogen (N)	[7727-37-9]	<0,0001	%		
	Nickel (Ni)	[7440-02-0]	0,55	%		
	Oxygen (O)	[7782-44-7]	~0,001	%		
	Phosphorus (P)	[7723-14-0]	0,006	%		
	Lead (Pb)	[7439-92-1]	<0,005	%		
	Sulfur (S)	[7704-34-9]	<0,001	%		
	Antimony (Sb)	[7440-36-0]	<0,005	%		
	Silicon (Si)	[7440-21-3]	0,033	%		
	Tin (Sn)	[7440-31-5]	0,01	%		
	Zinc (Zn)	[7440-66-6]	0,014	%		