

Copper base, Nickel Alloy, CRM, 706, Chips 100 g

Art. ID IARM-84C-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	~0,002	%		
	Carbon (C)	[7440-44-0]	~0,003	%		
	Cobalt (Co)	[7440-48-4]	~0,006	%		
	Copper (Cu)	[7440-50-8]	89,4	%		
	Iron (Fe)	[7439-89-6]	0,28	%		
	Manganese (Mn)	[7439-96-5]	~0,04	%		
	Nickel (Ni)	[7440-02-0]	9,8	%		
	Phosphorus (P)	[7723-14-0]	~0,06	%		
	Lead (Pb)	[7439-92-1]	0,047	%		
	Sulfur (S)	[7704-34-9]	~0,006	%		
	Antimony (Sb)	[7440-36-0]	~0,0012	%		
	Silicon (Si)	[7440-21-3]	~0,004	%		
	Tin (Sn)	[7440-31-5]	0,13	%		
	Zinc (Zn)	[7440-66-6]	0,13	%		