

**Copper base, Bronze - Leaded, CRM, 937, Disc 31 mm x 3 mm**

Art. ID IARM-92C-X  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	~0,05	%		
	Aluminium (Al)	[7429-90-5]	0,0013	%		
	Arsenic (As)	[7440-38-2]	~0,005	%		
	Beryllium (Be)	[7440-41-7]	<0,001	%		
	Bismuth (Bi)	[7440-69-9]	~0,011	%		
	Carbon (C)	[7440-44-0]	~0,002	%		
	Cobalt (Co)	[7440-48-4]	~0,0007	%		
	Chromium (Cr)	[7440-47-3]	~0,0007	%		
	Copper (Cu)	[7440-50-8]	80,35	%		
	Iron (Fe)	[7439-89-6]	~0,008	%		
	Manganese (Mn)	[7439-96-5]	~0,0016	%		
	Nickel (Ni)	[7440-02-0]	0,17	%		
	Oxygen (O)	[7782-44-7]	~0,004	%		
	Phosphorus (P)	[7723-14-0]	0,073	%		
	Lead (Pb)	[7439-92-1]	9,42	%		
	Sulfur (S)	[7704-34-9]	0,026	%		
	Antimony (Sb)	[7440-36-0]	0,078	%		
	Selenium (Se)	[7782-49-2]	~0,001	%		
	Silicon (Si)	[7440-21-3]	~0,0019	%		
	Tin (Sn)	[7440-31-5]	9,65	%		
	Zinc (Zn)	[7440-66-6]	0,146	%		