

**Iron base, Steel - Stainless, CRM, 410, Chips 100 g**

Art. ID IARM-FE410-18-C  
Unit 100 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,003	%		
	Carbon (C)	[7440-44-0]	0,132	%		
	Cobalt (Co)	[7440-48-4]	0,012	%		
	Chromium (Cr)	[7440-47-3]	12,2	%		
	Copper (Cu)	[7440-50-8]	0,046	%		
	Manganese (Mn)	[7439-96-5]	0,50	%		
	Molybdenum (Mo)	[7439-98-7]	0,146	%		
	Nitrogen (N)	[7727-37-9]	0,046	%		
	Niobium (Nb)	[7440-03-1]	0,0021	%		
	Nickel (Ni)	[7440-02-0]	0,28	%		
	Oxygen (O)	[7782-44-7]	0,009	%		
	Phosphorus (P)	[7723-14-0]	0,017	%		
	Sulfur (S)	[7704-34-9]	0,0014	%		
	Silicon (Si)	[7440-21-3]	0,29	%		
	Vanadium (V)	[7440-62-2]	0,065	%		
	Tungsten (W)	[7440-33-7]	0,008	%		
	Aluminium (Al)	[7429-90-5]	~0,003	%		
	Carbon (C)	[7440-44-0]	0,132	%		
	Cobalt (Co)	[7440-48-4]	0,012	%		
	Chromium (Cr)	[7440-47-3]	12,2	%		
	Copper (Cu)	[7440-50-8]	0,046	%		
	Manganese (Mn)	[7439-96-5]	0,5	%		
	Molybdenum (Mo)	[7439-98-7]	0,146	%		
	Nitrogen (N)	[7727-37-9]	0,046	%		
	Niobium (Nb)	[7440-03-1]	0,0021	%		
	Nickel (Ni)	[7440-02-0]	0,28	%		
	Oxygen (O)	[7782-44-7]	~0,009	%		
	Phosphorus (P)	[7723-14-0]	0,017	%		
	Sulfur (S)	[7704-34-9]	0,0014	%		
	Silicon (Si)	[7440-21-3]	0,29	%		
	Vanadium (V)	[7440-62-2]	0,065	%		

Tungsten (W)

[7440-33-7]

~0,008

%