

Nickel Alloy, Nickel Alloy 25-6Mo, UNS N08926, ~ 31 x 18 mm

Art. ID IARM-NI256MO-18-G
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,17 | % | | |
| | Arsenic (As) | [7440-38-2] | <0,0020 | % | | |
| | Boron (B) | [7440-42-8] | 0,0009 | % | | |
| | Carbon (C) | [7440-44-0] | 0,0182 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,02 | % | | |
| | Chromium (Cr) | [7440-47-3] | 20,3 | % | | |
| | Copper (Cu) | [7440-50-8] | 1,98 | % | | |
| | Iron (Fe) | [7439-89-6] | 44,0 | % | | |
| | Magnesium (Mg) | [7439-95-4] | ~0,0092 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,83 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 6,63 | % | | |
| | Nitrogen (N) | [7727-37-9] | 0,022 | % | | |
| | Niobium (Nb) | [7440-03-1] | 0,010 | % | | |
| | Nickel (Ni) | [7440-02-0] | 25,5 | % | | |
| | Oxygen (O) | [7782-44-7] | ~0,0030 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,004 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,0011 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,41 | % | | |
| | Tin (Sn) | [7440-31-5] | ~0,0020 | % | | |
| | Tantalum (Ta) | [7440-25-7] | <0,0300 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,006 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,009 | % | | |
| | Tungsten (W) | [7440-33-7] | 0,010 | % | | |
| | Zirconium (Zr) | [7440-67-7] | ~0,0030 | % | | |
| | Calcium (Ca) | [7440-70-2] | 0,0014 | % | | |
| | Aluminium (Al) | [7429-90-5] | 0,17 | % | | |
| | Boron (B) | [7440-42-8] | 0,0009 | % | | |
| | Carbon (C) | [7440-44-0] | 0,0182 | % | | |
| | Calcium (Ca) | [7440-70-2] | 0,0014 | % | | |
| | Cobalt (Co) | [7440-48-4] | ~0,02 | % | | |
| | Chromium (Cr) | [7440-47-3] | 20,3 | % | | |

| | | | |
|-----------------|-------------|--------|---|
| Copper (Cu) | [7440-50-8] | 1,98 | % |
| Iron (Fe) | [7439-89-6] | 44 | % |
| Manganese (Mn) | [7439-96-5] | 0,83 | % |
| Molybdenum (Mo) | [7439-98-7] | 6,63 | % |
| Nitrogen (N) | [7727-37-9] | 0,022 | % |
| Niobium (Nb) | [7440-03-1] | ~0,010 | % |
| Nickel (Ni) | [7440-02-0] | 25,5 | % |
| Phosphorus (P) | [7723-14-0] | ~0,004 | % |
| Sulfur (S) | [7704-34-9] | 0,0011 | % |
| Silicon (Si) | [7440-21-3] | 0,41 | % |
| Titanium (Ti) | [7440-32-6] | ~0,006 | % |
| Vanadium (V) | [7440-62-2] | ~0,009 | % |
| Tungsten (W) | [7440-33-7] | ~0,010 | % |