

17 components in methanol, CRM, ISO 17034 and ISO 17025

Art. ID LM24-CUS-70315
Unit 20 x 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methyl-tert.butylether (MTBE)	[1634-04-4]	40	mg/L		
	tert-Amyl-methyl ether	[994-05-8]	40	mg/L		
	2-Methylpentane	[107-83-5]	40	mg/L		
	2-Ethyltoluene	[611-14-3]	40	mg/L		
	3-Ethyltoluene	[620-14-4]	40	mg/L		
	1,2,3-Trimethylbenzene	[526-73-8]	40	mg/L		
	1,3,5-Trichlorobenzene	[108-70-3]	40	mg/L		
	Hexachloroethane	[67-72-1]	40	mg/L		
	Allyl chloride	[107-05-1]	40	mg/L		
	3-Chlorotoluene	[108-41-8]	40	mg/L		
	alpha-Chlorotoluene	[100-44-7]	40	mg/L		
	Bromotrichloromethane	[75-62-7]	40	mg/L		
	1-Bromo-2-chloroethane	[107-04-0]	40	mg/L		
	1,1,2,2-Tetrabromoethan e	[79-27-6]	40	mg/L		
	4-Ethyltoluene	[622-96-8]	40	mg/L		
	Ethyl-tert-butylether (ETBE)	[637-92-3]	40	mg/L		
	n-Hexane	[110-54-3]	40	mg/L		