

Custom Solution - 25 components in HNO₃ (7%), CRM, ISO 17034 and ISO 17025

Art. ID LM24-CUS-70322
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 12 months. (6 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cadmium (Cd)	[7440-43-9]	1	mg/L		
	Silver (Ag)	[7440-22-4]	5	mg/L		
	Barium (Ba)	[7440-39-3]	5	mg/L		
	Bismuth (Bi)	[7440-69-9]	5	mg/L		
	Cobalt (Co)	[7440-48-4]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	5	mg/L		
	Copper (Cu)	[7440-50-8]	5	mg/L		
	Gallium (Ga)	[7440-55-3]	5	mg/L		
	Lithium (Li)	[7439-93-2]	5	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Strontium (Sr)	[7440-24-6]	5	mg/L		
	Vanadium (V)	[7440-62-2]	5	mg/L		
	Boron (B)	[7440-42-8]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		
	Iron (Fe)	[7439-89-6]	10	mg/L		
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	10	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Phosphorus (P)	[7723-14-0]	20	mg/L		
	Magnesium (Mg)	[7439-95-4]	50	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		