

Custom Solution - 6 components of in Dichloromethane

Art. ID LM24-CUS-70324
Unit 5 x 1 mL
Deliverydetails De Minimis UN1593 III E1

Description

Conservation (closed) 12 months.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Naphthalene D8	[1146-65-2]	100	µg/mL		
	Fluorene D10	[81103-79-9]	100	µg/mL		
	Benzo(g,h,i)perylene D1 2	[93951-66-7]	100	µg/mL		
	Pyrene D10	[1718-52-1]	100	µg/mL		
	Benzo(a)pyrene D12	[63466-71-7]	100	µg/mL		
	Phenanthrene D10	[1517-22-2]	100	µg/mL		