

**Levothyroxine Peak Identification Mixture (1 mg/mL) (Levothyroxine sodium spiked with liothyronine, triiodothyroacetic acid and tetraiodothyroacetic acid in methanol) (COLD SHIPMENT REQUIRED)**

Art. ID                    USP-1365043  
Unit                        3 x 5 mL  
Deliverydetails         Dry ice shipment UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
(L-Tyrosine, O-(4-hydroxy-3,5-diiodophenyl)-3,5-diiodo-, monosodium salt)	Levothyroxine Sodium	[55-03-8]				
(O-(4-hydroxy-3-iodophenyl)-3,5-diiodo-) L-Tyrosine)	Liothyronine	[6893-02-3]				
([4-(4-Hydroxy-3-iodophenoxy)-3,5-diiodophenyl]acetic acid)	Triiodothyroacetic Acid	[51-24-1]				
(2-(4-(4-Hydroxy-3,5-diodophenoxy)-3,5-diiodophenyl)acetic acid)	Tetraiodothyroacetic Acid	[67-30-1]				