

## **Custom Solution - 8 components in HNO3 5%, CRM, ISO 17034 and ISO 17025**

Art. ID LM24-CUS-70338  
Unit 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### Description

Conservation 24 months. (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	400	mg/L		
	Arsenic (As)	[7440-38-2]	200	mg/L		
	Boron (B)	[7440-42-8]	800	mg/L		
	Chromium (Cr)	[7440-47-3]	800	mg/L		
	Iron (Fe)	[7439-89-6]	800	mg/L		
	Manganese (Mn)	[7439-96-5]	800	mg/L		
	Nickel (Ni)	[7440-02-0]	800	mg/L		
	Zinc (Zn)	[7440-66-6]	200	mg/L		