

Custom Solution - 2 components in HNO3 5%, 1% HF, CRM, ISO 17034 and ISO 17025

Art. ID LM24-CUS-70339
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Conservation 24 months. (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	8000	mg/L		
	Tin (Sn)	[7440-31-5]	4000	mg/L		