

Copper Alloy / CDA 642 / UNS C64200, CRM, Chips 100 g

Art. ID IARM-Cu642-18-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	6,5 ± 0,1	%		
	Copper (Cu)	[7440-50-8]	90,0 ± 1,0	%		
	Iron (Fe)	[7439-89-6]	0,039 ± 0,002	%		
	Manganese (Mn)	[7439-96-5]	0,0024 ± 0,0005	%		
	Nickel (Ni)	[7440-02-0]	0,014 ± 0,003	%		
	Phosphorus (P)	[7723-14-0]	0,009 ± 0,004	%		
	Lead (Pb)	[7439-92-1]	0,019 ± 0,002	%		
	Silicon (Si)	[7440-21-3]	1,96 ± 0,09	%		
	Tin (Sn)	[7440-31-5]	0,019 ± 0,002	%		
	Zinc (Zn)	[7440-66-6]	1,19 ± 0,05	%		