

Phosphatidylcholine, PC (soybean, hydrogenated), Purity: >98%

Art. ID LA-37-0108-0
Unit on request
Deliverydetails No Dangerous Good /not restricted

Description

Synonyms: PC (soybean, hydrogenated) / L-?-phosphatidylcholine (soybean, hydrogenated)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylcholine, PC (soybean, hydrogenated)	[97281-48-6]				