

**Phosphatidylinositol (liver, bovine), Purity: >98%, solution in Chloroform, concentration: Varies**

Art. ID LA-37-0134-0  
Unit on request  
Deliverydetails No Dangerous Good UN1888 III E1

Description

Synonyms: PI (liver, bovine) / L-?-phosphatidylinositol (liver, bovine)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylinositol (liver, bovine)	[383907-33-3]	Varies			