

Custom Solution - 3 components 1000 µg/mL each in Methanol

Art. ID LM24-CUS-70377
Unit 5 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ethylene Glycol	[107-21-1]	1000	µg/mL		
	Diethylene glycol	[111-46-6]	1000	µg/mL		
	Triethylene glycol	[112-27-6]	1000	µg/mL		