

Mineral Oil - Mix 4 in n-Hexane

Art. ID LM24-CUS-10093
Unit 1 x 1 mL
Deliverydetails De Minimis UN1208 II E2

Description

Conservation 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mineral oil Type A	[n/a]	5000	µg/mL		
	Mineral oil Type B	[n/a]	5000	µg/mL		
	n-Tetracontane (C40)	[4181-95-7]	2	µg/mL		
	n-Decane (C10)	[124-18-5]	1,46	µg/mL		