

## **Custom Solution - 5 components in Methanol / Dichloromethane (1/1)**

Art. ID LM24-CUS-70393  
Unit 3 x 1 mL  
Deliverydetails De Minimis UN1922 (6.1) II E2

### Description

Conservation 12 months (originally closed bottle)

| Text/Information | Analyte/Parameter  | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|--------------------|------------|---------------------|-------|--------|--------|
|                  | Dibenzo(a,e)pyrene | [192-65-4] | 100                 | µg/mL |        |        |
|                  | Dibenzo(a,h)pyrene | [189-64-0] | 100                 | µg/mL |        |        |
|                  | Dibenzo(a,i)pyrene | [189-55-9] | 100                 | µg/mL |        |        |
|                  | Dibenzo(a,l)pyrene | [191-30-0] | 100                 | µg/mL |        |        |
|                  | Perylene           | [198-55-0] | 100                 | µg/mL |        |        |