

Saccharin-d4 10 µg/mL in Methanol

Art. ID LM24-N-18000-1690D-10ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|----------------|---------------------|-------|--------|--------|
| | Saccharin-d4 | [1189466-17-8] | 10 | µg/mL | | |