

**n-Tetracosane-d50 (C24-d50) 100 µg/mL in 2,2,4-Trimethylpentane (Isooctane)**

Art. ID LM24-N-14700-0731D-100IO5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Tetracosane (D50,98%)	[16416-32-3]	100	µg/mL		