

### **Gondoic acid 1000 µg/mL in Acetone**

Art. ID LM24-N-20000-0126-1000AC5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

| Text/Information | Analyte/Parameter  | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|--|-------------|---------------------|-------|--------|--------|
|                  | (Z)-11-Eicosenoic acid<br>(Gondoic acid) (C20:1 n<br>-9) | [5561-99-9] | 1000                | µg/mL |        |        |